



enzyQUEST

Engineered M-MLV Reverse Transcriptase Basic Kit, 50 rxns Cat No: RN012S

Description:

Engineered M-MLV Reverse Transcriptase Basic Kit include a rationally designed M-MLV RT that exhibits increased thermal stability and reduced RNase H activity (RNase H⁻) than the wild type M-MLV RT.

- Engineered M-MLV RT remains active up to 50°C increasing the length and yield of cDNA compared to WT version.
- Able to efficiently detect and produce cDNA in reactions that have less than 100 pg of RNA template.
- Able to amplify at least 8kb in length

Source:

Gene coding the respective variant is heterologously expressed in *E.coli*.

Unit Definition:

One unit is defined as the amount of enzyme needed to catalyze the incorporation of 1 nmol dTTP into acid-insoluble material in 10 min under 37°C with poly(rA) and oligo(dT).

Reagents supplied:

| Cat number: RN012S | | |
|--------------------|--|-----------|
| Cat No | Description | Tubes No. |
| RN001S | Engineered MMLV Reverse transcriptase, 200u/μl, 0.05ml (200 units/μl in 20 mM Tris-HCl pH 7.4, 100 mM NaCl, 0.1 mM EDTA, 50 % glycerol, 0.01% IGEPAL, 1mM DTT) | 1 |
| BR055-1 | 5x Reverse Transcriptase buffer, 1ml (250 mM Tris-HCl (pH 8.3), 375 mM KCl, 15 mM MgCl ₂) | 1 |
| BR056-0.5 | 100mM DTT, 0.5ml | 1 |
| BR028-0.15 | dNTPs mix 10mM each, 0.15ml | 1 |
| BR033-1 | Nuclease Free Water, 1 ml | 1 |



Reaction guidelines and tips:

Quantities per Reaction (20 μ L Final Reaction Volume):

| | |
|---|--|
| 1. Preparation of Mix A | <i>1 ng -2 μg of total RNA*</i> <i>2 μL oligod(T)₁₈ from 50 μM stock §, §§</i> <i>1 μL dNTPs mix from 10 mM each stock</i> <i>X μL of Nuclease Free Water to 10 μL</i> |
| 2. Incubate Mix A at 65° C for 5 min, then transfer for 1 min on ice (meanwhile prepare Mix B) | |
| 3. Preparation of Mix B | <i>2,5 μL of Nuclease Free Water</i> <i>4 μL from 5x RT Buffer</i> <i>2 μL DTT from 100 mM stock</i> <i>0,5 μL of RNase inhibitor**</i> <i>1 μL of Reverse Transcriptase enzyme from stock of 200 U/μL</i> |
| 4. Following Step 2, add 10 μL of Mix B to the reaction | |
| 5. Incubate at 37° C to 55° C for 15-60 minutes *** | |
| 6. Heat inactivate the Reverse Transcriptase at 65° C for 20 min | |
| 7. Cool down the cDNA to room temperature and use 0,1 - 2 μL in PCR reaction | |

* $\leq 0,1$ ng may suffice for detecting certain RNA populations that are abundant like e.g. Rubisco transcript in plants. Use 10 times less RNA quantity for coding RNAs enriched substrates.

** The final concentration of the RNase Inhibitor is 1 U/ μ L. Since RT enzyme and reagents have been checked for the presence of RNases and DNases, this step may be omitted if the sample is of high RNA quality.

*** Routine temp for a standard reaction is 42° C. For most RNA cases tested complete cDNA production was verified within **10-15 minutes** independently of the incubation temperature. [Optimize depending on your target and provide us with your feedback to get special offers and rewards!](#)

§ In case of random hexamers, add 2 μ L from 50 μ M stock per 20 μ L RT reactions. Following “step2” incubate for 5 minutes at 25° C or room temperature before incubating the reaction to 42° C.

§§ In case of gene specific primer, add 2 μ L from 2 μ M primer stock. Incubate at 50° C for cDNA production. Optimize to your case, depending on primer's T_m. Incorporate 5-10 degrees below the T_m of your oligo. Do not incubate at temperatures higher than 55° C.



enzyQUEST

Functional Quality Control:

Two step cDNA synthesis reaction in order to amplify an 831-bp region using 1ng Hela RNA as substrate. The resulting PCR product is visualized as a single band on a midori-stained agarose gel.

Other Quality Controls

Tested extensively for the absence of DNases and RNases and *E.coli* DNA contamination.

Shipping

Shipped on blue ice or dry ice

Storage conditions

Store at $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$

Shelf life:

24 months post production

Version 3.0/ Jan-26